




THUASNE

BIFLEX Self Adjust

fr	Dispositif compressif ajustable pour mollet, cuisse et bras	8
en	Adjustable compression device for calf, thigh and arm	8
de	Verstellbare Kompressionsbandage für Wade, Oberschenkel und Arm	9
nl	Verstelbaar compressiehulpmiddel voor kuit, dij en arm	10
da	Justerbar kompressionsanordning til læg, lår og arm	10
sv	Justerbar kompressionsanordning för vad, lår och arm	11
no	Justerbar kompresjonsinnretning for legg, lår og arm	12

Mollet - Calf - Wade - Kuit - Læg - Vad - Legg

	ℓ	cB	cC	cD	N	L
fr	Hauteur cheville-genou	Circonférence de cheville (au plus fin)	Circonférence de mollet (au plus fort)	Circonférence sous le genou	Normal	Long
en	Ankle-knee height	Ankle circumference (narrowest point)	Calf circumference (widest point)	Circumference below the knee	Normal	Long
de	Knöchel-Knie-Höhe	Knöchelumfang (an der schmalsten Stelle)	Wadenumfang (an der stärksten Stelle)	Umfang unterhalb des Knies	Normal	Lang
nl	Enkel-knie hoogte	Omtrek van de enkel (kleinste omtrek)	Omtrek van de kuit (grootste omtrek)	Omtrek onder de knie	Normaal	Lang
da	Ankel-knæ højde	Anklens omkreds (det tyndeste sted)	Læggens omkreds (det kraftigste sted)	Omkreds under knæ	Normal	Lang
sv	Höjd ankel-knä	Omkrets ankel (på det smalaste stället)	Omkrets vad (på det bredaste stället)	Omkrets under knä	Normal	Lång
no	Ankel-kne høyde	Anklens omkrets (på det tynneste)	Leggens omkrets (på det tykkeste)	Omkrets under kneet	Normal	Lang

	cB	cC	cD
1	16 - 24 cm	24 - 35 cm	25 - 36 cm
	6.5 - 9.5"	9.5 - 14"	10 - 14"
2	20 - 28 cm	32 - 43 cm	32 - 43 cm
	8 - 11"	12.5 - 17"	12.5 - 17"
3	24 - 32 cm	40 - 51 cm	39 - 50 cm
	9.5 - 12.5"	15.5 - 20"	15.5 - 19.5"
4	28 - 36 cm	48 - 59 cm	46 - 57 cm
	11 - 14"	19 - 23"	18 - 22.5"

N	L
ℓ < 36 cm / 14"	ℓ > 36 cm / 14"

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- en**
1. Before putting it on, ensure that the hook-and-loop fasteners are rolled up. ①
 2. Slip the arm into the sleeve, Thuasne logo upwards, sleeve on the inside of the arm. ②
 3. Place the bottom ④ strap: hold the inner part against the skin, tighten the other part over it and fasten with the hook-and-loop fastener (do not fasten to the sleeve). ③
 4. Unroll the other straps from bottom to top ⑥, ⑦, ⑧, ⑨.
 5. If necessary, wrap the removable strap over the elbow. ⑩
The hook-and-loop fastener must not be closed over the sleeve.
 6. To remove it, undo the straps from top to bottom and roll the hook-and-loop fasteners up again and pull the product off. ⑪

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Conservé cette notice.

en

ADJUSTABLE COMPRESSION DEVICE FOR CALF, THIGH AND ARM

Description/Destination

The device comprises a number of adjustable and repositionable overlapping straps which are closed with hook-and-loop fasteners.

Different models are available to fit various parts of the body: foot, calf, knee, thigh, arm and hand.

The calf and thigh models have a tightening system used to control and reproduce the desired level of pressure: low (minimum 20 mmHg), medium (minimum 30 mmHg) and high (minimum 40 mmHg).

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Textile components: polyamide - elastane - polyurethane.

Tightening system: polyamide - polyurethane.

Properties/Mode of action

The devices provide compression by applying pressure to the limb. They promote venous and lymphatic circulation, help prevent the formation of edema and maintain volume.

Inelastic short-stretch material (average pressure between 20 and 50 mmHg).

Indications

Management of edema (venous, lymphatic and caused by injury).

Chronic venous insufficiency.

Symptomatic treatment of lipedema.

Contraindications

Do not use in case of major dermatological conditions on an affected limb.

Do not use in the event of known allergy to any of the components.

Do not use in case of lower limb peripheral arterial disease (PAD) with ankle-brachial index (ABI) < 0.6.

Do not use in the event of decompensated heart failure.

Do not use in the event of advanced diabetic microangiopathy (for a compression > 30 mmHg).

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in the event of phlegmasia cerulea dolens (painful blue phlebitis with arterial compression) of an affected limb.

Do not use in case of extra-anatomical bypass of the affected limb.

Do not use in the event of septic thrombosis.

Do not use on the upper limbs in the following cases:

- Do not use in case of brachial plexus condition;
- Do not use in case of vasculitis of the extremities.

Precautions

Strictly comply with your healthcare professional's prescription and recommendations for use.

If the device seems to be too tight at the recommended tightness, remove it and see a healthcare professional.

The risk/benefit ratio and the adequate level of pressure must be regularly re-assessed by a healthcare professional in the event of:

- PAD with ABPI between 0.6 and 0.9
- Advanced peripheral neuropathy
- Weeping or eczematous dermatitis.

For hygiene, security and performance reasons, do not re-use the product for another patient.

In the event of discomfort, significant restriction, pain, variation in limb volume, deterioration of skin condition, infection, unusual sensations, change in the colour of the extremities, or changes in the performance of the product, remove the product and seek the advice of a healthcare professional.

Examine the shape of the leg: if there are any deformed or bony areas, use the protective foam bandages or Varico pads to make the shape of the limb uniform.

In the event of a venous ulcer, a Doppler examination is necessary before using the device in order to ensure there is no associated severe arterial disease.

In the event of lipedema, compression can be used to help relieve the symptoms. Seek advice from a healthcare professional for compression and other solutions for managing lipedema.

For some pathologies (or situations) such as venous thrombosis, the device must be used in combination with an anticoagulant treatment; consult a healthcare professional for advice.

Do not use products on the skin (creams, ointments, oils, gels, patches...) before applying the device, it could damage the product.

Do not apply the product directly to an open wound without a dressing.

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity or cutaneous dryness.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

- **Putting the calf model on:** pages 2/3

- **Putting the thigh model on:** pages 4/5

- **Putting the arm model on:** pages 6/7

When putting it on and wearing it, ensure no large folds appear.

If necessary, the segment pressure can be adjusted by adjusting the loops from bottom to top.

Care/Maintenance

Compression devices are designed for daily use over a 6-month period (i.e. 90 washes). After this time, garments should be replaced to ensure the treatment is effective and technical performance is not affected.

Close the hook-and-loop fasteners before washing. Machine washable at 30 °C (delicate programme). If possible use a washing net. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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©Thuasne - 2043201 (2022-12)
THUASNE UK Ltd