# THUASNE

## Genuimmo 20°

- fr** Attelle d'immobilisation du genou à 20°
- en** Knee immobilisation splint at 20°
- de** Schiene zur Ruhigstellung des Knies in einem 20°-Winkel
- nl** Immobilisatiebrace voor de knie op 20°
- it** Tutore di immobilizzazione per il ginocchio a 20°
- es** Férula de inmovilización de la rodilla a 20°
- pt** Tala para imobilização do joelho em 20°
- da** Immobiliserende knæskinne ved 20°
- fi** Immobilisoiva polvituki, taivutus 20°
- sv** Immobiliseringssskena för knät 20°
- el** Νάρθηκας ακινητοποίησης γονάτου σε κάμψη 20°
- cs** Imobilizační ortéza kolene 20°
- pl** Orteza do unieruchomienia kolana pod kątem 20°
- lv** Ceļgala imobilizācijas sāite 20°
- lt** Kelio imobilizavimo 20° padėtyje įtvaras
- et** Põlve fikseerimise lahas 20° juures
- sl** Opornica za imobilizacijo kolena s kotom 20°
- sk** Imobilizačná kolenná ortéza pri 20° uhle
- hu** Oldalsínes térdrögzítő 20°-os szög
- bg** Шина за имобилизация на коляно на 20°
- ro** Atelă pentru imobilizarea genunchiului la 20°
- ru** Шина для имобилизации коленного сустава под углом 20°
- hr** Udlaga za imobilizaciju koljena, 20°
- zh** 20° 膝关节固定夹板
- ar** جبيرة لتثبيت الركبة على زاوية 20 درجة

		< 175 cm	> 175 cm
	<b>S</b>	35 - 45 cm	x
	<b>M</b>	45 - 52 cm	
	<b>L</b>	53 - 62 cm	
	<b>XL</b>	63 - 75 cm	60 cm




**THUASNE**  
120, rue Marius AUFAN  
92300 Levallois-Perret - France

© Thuasne - 6165201 (2020-12) © Studio Caterin

[www.thuasne.com](http://www.thuasne.com)

DM210444

		< 175 cm	> 175 cm
	<b>S</b>	35 - 45 cm	x
	<b>M</b>	45 - 52 cm	
	<b>L</b>	53 - 62 cm	
	<b>XL</b>	63 - 75 cm	
		50 cm	60 cm

<b>fr</b>	Taille patient	Circonférence cuisse	Hauteur	Immobilisation
<b>en</b>	Patient height	Thigh circumference	Height	Immobilisation
<b>de</b>	Größe des Patienten	Oberschenkelumfang	Höhe	Ruhigstellung
<b>nl</b>	Maat van de patiënt	Dijomtrek	Hoogte	Immobilisatie
<b>it</b>	Altezza del paziente	Circonferenza coscia	Lunghezza	Immobilizzazione
<b>es</b>	Tamaño del paciente	Circunferencia del muslo	Altura	Inmovilización
<b>pt</b>	Tamanho do paciente	Circunferência da coxa	Altura	Imobilização
<b>da</b>	Patientens højde	Lårets omkreds	Højde	Støtte
<b>fi</b>	Potilaan koko	Reiden ympärysmitta	Korkeus	Immobilisoi
<b>sv</b>	Patientens storlek	Omkrets runt lår	Höjd	Immobilisering
<b>el</b>	Μέγεθος ασθενούς	Περιφέρεια μηρού	Ύψος	Ακινητοποίηση
<b>cs</b>	Výška pacienta	Obvod stehna	Délka	Znehybnění
<b>pl</b>	Rozmiar pacjenta	Obwód uda	Wysokość	Unieruchomienie
<b>lv</b>	Pacienta augums	Augšstilba apkārtmērs	Augstums	Imobilizācija
<b>lt</b>	Paciento dydis	Šlaunies apimtis	Aukštis	Imobilizavimas
<b>et</b>	Patsiendi suurus	Reie ümbermõõt	Kõrgus	Liikumatuks muutmine
<b>sl</b>	Velikost pacienta	Obseg stegna	Višina	Imobilizacija
<b>sk</b>	Veľkosť pacienta	Obvod stehna	Výška	Imobilizácia
<b>hu</b>	A beteg mérete	Comb körmérete	Magasság	Rögzítés
<b>bg</b>	Размер на пациента	Бедрена обиколка	Височина	Обездвиждане
<b>ro</b>	Talia pacientului	Circumferința coapsei	Înălțime	Imobilizare
<b>ru</b>	Размер пациента	Окружность бедра	Высота	Иммобилизация
<b>hr</b>	Veličina pacijenta	Promjer bedra	Visina	Imobilizacija
<b>zh</b>	患者身高	大腿圍	高度	固定
<b>ar</b>	مقاس المريض	محيط الفخذ	الارتفاع	التثبيت

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Conservser cette notice.

**en**

## KNEE IMMOBILISATION SPLINTS AT 0° AND 20°

**Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

The device is composed of :

•3 removable, radiolucent, rigid stays.

•4 self-fastening closure and support straps.

Option of opening the side of the splint while still keeping the knee immobilised (see application).

Knee splint with fine fabric and well-ventilated thanks to large opening.

No seams against the skin.

Bilateral product.

Available in 4 sizes for the short version and 3 sizes for the long version.

**Composition**

Texile components: PVC – Polyamide – Cotton – Polyurethane – Polyester – Natural latex – Elastane.

Rigid components: Aluminium – Silicone – Acetal.

**Properties/Mode of action**

Knee joint stabilisation via the rigid stays.

Total immobilisation at 0° or 20°.

Adjustable closure and support via self-fastening straps.

**Indications**

Post-operative immobilisation.

Post-traumatic immobilisation.

Knee pain management.

**Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprohylaxis.

Do not use in the event of unstable fractures, fractures of the lower part of the tibia or fibula, or fractures of the upper part of the femur.

**Precautions**

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Prolonged immobilisation can sometimes cause muscle weakness.

Do not wear the product when driving a vehicle.

Contains natural latex.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

**Undesirable side-effect**

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

**Instructions for use/Application**

Choose the splint height, then the size.

Measurements should be taken at the widest part of the thigh.

It is recommended that a healthcare professional supervises the first application.

It is recommended to wear the device directly on the skin, unless contraindicated.

•Open the splint and lay it flat open.

•Place the opened splint under the leg, the wider edge towards the top of the leg.

**N.B.**: The top of the splint must not impede flexion of the leg at the groin.

•Shaping the stays: if needed, it is possible to shape the stays in order to avoid potential pressure points or the splint slipping. Apply a slight pressure in the desired direction.

Repeat the process on the other stays if needed.

If subsequent shaping of the stays is needed, consult a healthcare professional.

•Tighten the straps in the order described.

Tightening of ❶ and ❷ by passing them under the lateral uprights of the splint.

Fastening of the panels ❶ and ❷.

Tightening of ❸ and ❹ by passing throught the buckles.

Once all the straps are fastened, readjust the tightening.

•Splint adaptability:

The straps located just above and below the knee are height-adjustable ❺, meaning that they can be adapted to fit the patient's shape and condition (presence of scar following surgery, wound, etc.).

The knee splint can be opened on the outside ❻, particularly in the event of surgery.

To do this, simply unfasten the hook located in the middle of the splint and pass it under the upright before fastening it at the back of the knee splint.

If needed, reduce the strap length by cutting the edges with scissors.

**Care**

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. Machine washable at 30 °C (delicate programme). Close the self-fastening tabs before washing. Remove stays before washing. Stays can be washed with a moist cloth.

Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

**Storage**

Store at room temperature, preferably in the original packaging.

**Disposal**

Dispose of in accordance with local regulations.

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